



budC Antibody

Product Code	CSB-PA166347XA01FRW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ZNN8
Immunogen	Recombinant Corynebacterium glutamicum (Brevibacterium saccharolyticum) budC protein
Raised In	Rabbit
Species Reactivity	Corynebacterium glutamicum (Brevibacterium saccharolyticum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Corynebacterium glutamicum (Brevibacterium saccharolyticum)
Lead Time	Made-to-order (14-16 weeks)
Target Names	budC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.