



Capsid Antibody

Product Code	CSB-PA361107XA01ECK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03614
Immunogen	Recombinant Enterobacteria phage fr (Bacteriophage fr) Capsid protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage fr (Bacteriophage fr)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage fr (Bacteriophage fr)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.