



bgaA Antibody

Product Code	CSB-PA527526XA01TNQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O54315
Immunogen	Recombinant Thermus sp. bgaA protein
Raised In	Rabbit
Species Reactivity	Thermus sp.
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Thermus sp.
Lead Time	Made-to-order (14-16 weeks)
Gene Names	bgaA