



bnzB Antibody

Product Code	CSB-PA405973XA01FRZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5W4F1
Immunogen	Recombinant Pseudomonas putida (strain ATCC 700007 / DSM 6899 / BCRC 17059 / F1) bnzB protein
Raised In	Rabbit
Species Reactivity	Pseudomonas putida (strain ATCC 700007 / DSM 6899 / BCRC 17059 / F1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas putida (strain ATCC 700007 / DSM 6899 / BCRC 17059 / F1)
Lead Time	Made-to-order (14-16 weeks)
Target Names	bnzB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.