



bioY Antibody

Product Code	CSB-PA326347XA01LRU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22819
Immunogen	Recombinant Lysinibacillus sphaericus (Bacillus sphaericus) bioY protein
Raised In	Rabbit
Species Reactivity	Lysinibacillus sphaericus (Bacillus sphaericus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lysinibacillus sphaericus (Bacillus sphaericus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	bioY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.