



Outer membrane lipoprotein B1c Antibody

Product Code	CSB-PA669354XA01DTO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q46036
Immunogen	Recombinant Citrobacter freundii Outer membrane lipoprotein B1c protein
Raised In	Rabbit
Species Reactivity	Citrobacter freundii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Citrobacter freundii
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.