



caf1R Antibody

Product Code	CSB-PA340829XA01YAS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26950
Immunogen	Recombinant Yersinia pestis caf1R protein
Raised In	Rabbit
Species Reactivity	Yersinia pestis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Yersinia pestis
Lead Time	Made-to-order (14-16 weeks)
Target Names	caf1R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.