



Alpha-conotoxin LvIA Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA926267XA01DWC |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | L8BU87 |
| Immunogen | Recombinant Conus lividus (Livid cone) Alpha-conotoxin LvIA protein (21-36aa) |
| Raised In | Rabbit |
| Species Reactivity | Conus lividus (Livid cone) |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Conus lividus (Livid cone) |
| Target Names | Alpha-conotoxin LvIA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |