



gp24 Antibody

Product Code	CSB-PA191539XA01EEA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D2K805
Immunogen	Recombinant Enterobacteria phage T6 (Bacteriophage T6) gp24 protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T6 (Bacteriophage T6)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage T6 (Bacteriophage T6)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gp24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.