



cas1 Antibody

Product Code	CSB-PA570647XA01FNQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G3ECR2
Immunogen	Recombinant Streptococcus thermophilus cas1 protein
Raised In	Rabbit
Species Reactivity	Streptococcus thermophilus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptococcus thermophilus
Lead Time	Made-to-order (14-16 weeks)
Target Names	cas1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.