



CP Antibody

Product Code	CSB-PA361096XA01OAM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03578
Immunogen	Recombinant Odontoglossum ringspot virus (isolate Japan) (ORSV) CP protein
Raised In	Rabbit
Species Reactivity	Odontoglossum ringspot virus (isolate Japan) (ORSV)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Odontoglossum ringspot virus (isolate Japan) (ORSV)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.