



CYP98A1 Antibody

Product Code	CSB-PA529860XA01FJS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O48956
Immunogen	Recombinant Sorghum bicolor (Sorghum) (Sorghum vulgare) CYP98A1 protein
Raised In	Rabbit
Species Reactivity	Sorghum bicolor (Sorghum) (Sorghum vulgare)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sorghum bicolor (Sorghum) (Sorghum vulgare)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CYP98A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.