




Capsid protein Antibody

Product Code	CSB-PA356100XA01FQZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03615
Immunogen	Recombinant Escherichia phage Qbeta (Bacteriophage Q-beta) CAPSD protein
Raised In	Rabbit
Species Reactivity	Escherichia phage Qbeta (Bacteriophage Q-beta)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage Qbeta (Bacteriophage Q-beta)

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: Capsid protein antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 18.3 kDa

Observed band size: 20 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.