



C Antibody

Product Code	CSB-PA355964XA01HEK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03147
Immunogen	Recombinant Hepatitis B virus genotype D subtype adw (isolate United Kingdom/adyw/1979) (HBV-D) C protein
Raised In	Rabbit
Species Reactivity	Hepatitis B virus genotype D subtype adw (isolate United Kingdom/adyw/1979) (HBV-D)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Hepatitis B virus genotype D subtype adw (isolate United Kingdom/adyw/1979) (HBV-D)
Lead Time	Made-to-order (14-16 weeks)
Target Names	C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.