



katB Antibody

Product Code	CSB-PA950908XA01PWD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46206
Immunogen	Recombinant Pseudomonas syringae pv. syringae katB protein
Raised In	Rabbit
Species Reactivity	Pseudomonas syringae pv. syringae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas syringae pv. syringae
Lead Time	Made-to-order (14-16 weeks)
Target Names	katB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.