




cdd Antibody

Product Code	CSB-PA004976XA01ENP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1IYC6
Immunogen	Recombinant Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks) cdd protein (1-294aa)
Raised In	Rabbit
Species Reactivity	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli (strain ATCC 8739 / DSM 1576 / Crooks)
Target Names	cdd

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes:cdd Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 38.5 kDa

Observed band size: 42 kDa

Exposure time: 10s.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.