



cgt Antibody

Product Code	CSB-PA333219XA01BAW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27036
Immunogen	Recombinant Bacillus ohbensis cgt protein
Raised In	Rabbit
Species Reactivity	Bacillus ohbensis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus ohbensis
Lead Time	Made-to-order (14-16 weeks)
Gene Names	cgt