




cheB Antibody

Product Code	CSB-PA362107XA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07330
Immunogen	Synthesized peptide derived from internal of Escherichia coli (strain K12) cheB protein (246-260aa)
Raised In	Rabbit
Species Reactivity	Escherichia coli (strain K12)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli (strain K12)
Target Names	cheB

Image



Dot blot analysis of peptides corresponding to the immunogen was probed with cheB Antibody at a dilution of 1:1000. Blocking/Dilution buffer: 4% nonfat dry milk/PBS. Secondary antibody: Goat polyclonal to rabbit IgG at 1/40000 dilution Exposure time: 10s.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.