



CHI Antibody

Product Code	CSB-PA677099XA01PWZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q43056
Immunogen	Recombinant Pueraria montana var. lobata (Kudzu vine) (Pueraria lobata) CHI protein
Raised In	Rabbit
Species Reactivity	Pueraria montana var. lobata (Kudzu vine) (Pueraria lobata)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pueraria montana var. lobata (Kudzu vine) (Pueraria lobata)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CHI
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.