



groEL Antibody

Product Code	CSB-PA331509XA01EAW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42382
Immunogen	Recombinant Ehrlichia chaffeensis groEL protein
Raised In	Rabbit
Species Reactivity	Ehrlichia chaffeensis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ehrlichia chaffeensis
Lead Time	Made-to-order (14-16 weeks)
Target Names	groEL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.