



chy1 Antibody

Product Code	CSB-PA330215XA01GBX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32963
Immunogen	Recombinant Gibberella baccata (Fusarium lateritium) chy1 protein
Raised In	Rabbit
Species Reactivity	Gibberella baccata (Fusarium lateritium)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Gibberella baccata (Fusarium lateritium)
Lead Time	Made-to-order (14-16 weeks)
Target Names	chy1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.