



CDO Antibody

Product Code	CSB-PA866190XA01SXP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9U8F1
Immunogen	Recombinant Schistosoma japonicum (Blood fluke) CDO protein
Raised In	Rabbit
Species Reactivity	Schistosoma japonicum (Blood fluke)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Schistosoma japonicum (Blood fluke)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CDO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.