



CU Antibody

Product Code	CSB-PA313268XA01OBS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06153
Immunogen	Recombinant Ophiostoma ulmi (Dutch elm disease fungus) CU protein
Raised In	Rabbit
Species Reactivity	Ophiostoma ulmi (Dutch elm disease fungus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Ophiostoma ulmi (Dutch elm disease fungus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CU
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.