



groEL Antibody

Product Code	CSB-PA651539XA01FSS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2NW94
Immunogen	Recombinant Sodalis glossinidius (strain morsitans) groEL protein
Raised In	Rabbit
Species Reactivity	Sodalis glossinidius (strain morsitans)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sodalis glossinidius (strain morsitans)
Lead Time	Made-to-order (14-16 weeks)
Target Names	groEL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.