



CP84 Antibody

Product Code	CSB-PA620532XA01FRV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41062
Immunogen	Recombinant Escherichia phage 186 (Bacteriophage 186) CP84 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage 186 (Bacteriophage 186)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage 186 (Bacteriophage 186)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CP84
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.