



Cycloisomaltooligosaccharide glucanotransferase Antibody

Product Code	CSB-PA311318XA01BAN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P94286
Immunogen	Recombinant Bacillus circulans Cycloisomaltooligosaccharide glucanotransferase protein
Raised In	Rabbit
Species Reactivity	Bacillus circulans
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus circulans
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.