



coaE Antibody

Product Code	CSB-PA349716XA01DWW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58101
Immunogen	Recombinant Corynebacterium ammoniagenes (Brevibacterium ammoniagenes) coaE protein
Raised In	Rabbit
Species Reactivity	Corynebacterium ammoniagenes (Brevibacterium ammoniagenes)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Corynebacterium ammoniagenes (Brevibacterium ammoniagenes)
Lead Time	Made-to-order (14-16 weeks)
Target Names	coaE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.