



cry1Ab Antibody

Product Code	CSB-PA363359XA01BCU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A372
Immunogen	Recombinant Bacillus thuringiensis subsp. aizawai cry1Ab protein
Raised In	Rabbit
Species Reactivity	Bacillus thuringiensis subsp. aizawai
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus thuringiensis subsp. aizawai
Lead Time	Made-to-order (14-16 weeks)
Target Names	cry1Ab
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.