



cyt1Aa Antibody

Product Code	CSB-PA963496XA01BCZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A382
Immunogen	Recombinant Bacillus thuringiensis subsp. israelensis cyt1Aa protein
Raised In	Rabbit
Species Reactivity	Bacillus thuringiensis subsp. israelensis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus thuringiensis subsp. israelensis
Lead Time	Made-to-order (14-16 weeks)
Target Names	cyt1Aa
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.