



cya Antibody

Product Code	CSB-PA994876XA01BQE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40136
Immunogen	Recombinant Bacillus anthracis cya protein
Raised In	Rabbit
Species Reactivity	Bacillus anthracis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus anthracis
Lead Time	Made-to-order (14-16 weeks)
Target Names	cya
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.