



crtY Antibody

Product Code	CSB-PA322008XA01ERA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21687
Immunogen	Recombinant Pantoea ananas (Erwinia uredovora) crtY protein
Raised In	Rabbit
Species Reactivity	Pantoea ananas (Erwinia uredovora)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pantoea ananas (Erwinia uredovora)
Lead Time	Made-to-order (14-16 weeks)
Target Names	crtY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.