



ctaD Antibody

Product Code	CSB-PA886346XA01PCD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00502
Immunogen	Recombinant Paracoccus versutus (Thiobacillus versutus) ctaD protein
Raised In	Rabbit
Species Reactivity	Paracoccus versutus (Thiobacillus versutus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Paracoccus versutus (Thiobacillus versutus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ctaD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.