



cpnB Antibody

Product Code	CSB-PA810126XA01FSW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8GAW0
Immunogen	Recombinant Comamonas sp. (strain NCIMB 9872) cpnB protein
Raised In	Rabbit
Species Reactivity	Comamonas sp. (strain NCIMB 9872)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Comamonas sp. (strain NCIMB 9872)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	cpnB