



dcyD Antibody

Product Code	CSB-PA483395XA01EEO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7M363
Immunogen	Recombinant Escherichia coli O8 (strain IA11) dcyD protein
Raised In	Rabbit
Species Reactivity	Escherichia coli O8 (strain IA11)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli O8 (strain IA11)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	dcyD