



desA Antibody

Product Code	CSB-PA684488XA01FTE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q54794
Immunogen	Recombinant Arthrospira platensis (Spirulina platensis) desA protein
Raised In	Rabbit
Species Reactivity	Arthrospira platensis (Spirulina platensis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arthrospira platensis (Spirulina platensis)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	desA