



Cysteine synthase Antibody

Product Code	CSB-PA669274XA01CKU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q43317
Immunogen	Recombinant Citrullus lanatus (Watermelon) (Citrullus vulgaris) Cysteine synthase protein
Raised In	Rabbit
Species Reactivity	Citrullus lanatus (Watermelon) (Citrullus vulgaris)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Citrullus lanatus (Watermelon) (Citrullus vulgaris)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.