



crtB Antibody

Product Code	CSB-PA325131XA01ERA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21683
Immunogen	Recombinant Pantoea ananas (Erwinia uredovora) crtB protein
Raised In	Rabbit
Species Reactivity	Pantoea ananas (Erwinia uredovora)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pantoea ananas (Erwinia uredovora)
Lead Time	Made-to-order (14-16 weeks)
Target Names	crtB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.