




CRTI Antibody

Product Code	CSB-PA324241XA01ERA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21685
Immunogen	Recombinant Pantoea ananas (Erwinia uredovora) CRTI protein
Species Reactivity	Pantoea ananas (Erwinia uredovora)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pantoea ananas (Erwinia uredovora)
Target Names	CRTI

Image



Western Blot

Positive WB detected in Recombinant protein
 All lanes: CRTI antibody at 1:2000
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 24.2 kDa
 Observed band size: 27 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.