




hupA Antibody

Product Code	CSB-PA364791XA01FQR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0ACF1
Immunogen	Recombinant Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC) hupA protein (1-90aa)
Raised In	Rabbit
Species Reactivity	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Target Names	hupA Antibody

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes:hupA Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 19 kDa

Observed band size: 19 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.