



ntpT Antibody

Product Code	CSB-PA522466XA01LMP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O07380
Immunogen	Recombinant Lactobacillus helveticus (Lactobacillus suntoryeus) ntpT protein
Raised In	Rabbit
Species Reactivity	Lactobacillus helveticus (Lactobacillus suntoryeus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Lactobacillus helveticus (Lactobacillus suntoryeus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ntpT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.