



dnaJ Antibody

Product Code	CSB-PA486133XA01EOP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7MNM2
Immunogen	Recombinant Escherichia coli O81 (strain ED1a) dnaJ protein
Raised In	Rabbit
Species Reactivity	Escherichia coli O81 (strain ED1a)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli O81 (strain ED1a)
Lead Time	Made-to-order (14-16 weeks)
Target Names	dnaJ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.