



DERF3 Antibody

Product Code	CSB-PA341808XA01DCO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49275
Immunogen	Recombinant Dermatophagoides farinae (American house dust mite) DERF3 protein
Raised In	Rabbit
Species Reactivity	Dermatophagoides farinae (American house dust mite)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Dermatophagoides farinae (American house dust mite)
Lead Time	Made-to-order (14-16 weeks)
Target Names	DERF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.