



dex Antibody

Product Code	CSB-PA328009XA01SMJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P39653
Immunogen	Recombinant Streptococcus downei (Streptococcus sobrinus) dex protein
Raised In	Rabbit
Species Reactivity	Streptococcus downei (Streptococcus sobrinus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptococcus downei (Streptococcus sobrinus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	dex
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.