



dinG Antibody

Product Code	CSB-PA335488XA01FTM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29741
Immunogen	Recombinant Photobacterium profundum (strain SS9) dinG protein
Raised In	Rabbit
Species Reactivity	Photobacterium profundum (strain SS9)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Photobacterium profundum (strain SS9)
Lead Time	Made-to-order (14-16 weeks)
Target Names	dinG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.