



dgdA Antibody

Product Code	CSB-PA322324XA01BXR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16932
Immunogen	Recombinant Burkholderia cepacia (Pseudomonas cepacia) dgdA protein
Raised In	Rabbit
Species Reactivity	Burkholderia cepacia (Pseudomonas cepacia)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Burkholderia cepacia (Pseudomonas cepacia)
Lead Time	Made-to-order (14-16 weeks)
Target Names	dgdA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.