



foIP1 Antibody

Product Code	CSB-PA366228XA01FQW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0X1
Immunogen	Recombinant Mycobacterium leprae (strain TN) foIP1 protein
Raised In	Rabbit
Species Reactivity	Mycobacterium leprae (strain TN)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mycobacterium leprae (strain TN)
Lead Time	Made-to-order (14-16 weeks)
Target Names	foIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.