



folP1 Antibody

Product Code	CSB-PA782580XA01FSG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WND1
Immunogen	Recombinant Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv) folP1 protein
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Lead Time	Made-to-order (14-16 weeks)
Target Names	folP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.