



dsbC Antibody

Product Code	CSB-PA331060XA01DKD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P39691
Immunogen	Recombinant Dickeya dadantii (strain 3937) (Erwinia chrysanthemi (strain 3937)) dsbC protein
Raised In	Rabbit
Species Reactivity	Dickeya dadantii (strain 3937) (Erwinia chrysanthemi (strain 3937))
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Dickeya dadantii (strain 3937) (Erwinia chrysanthemi (strain 3937))
Lead Time	Made-to-order (14-16 weeks)
Target Names	dsbC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.