



DUT Antibody

Product Code	CSB-PA525044XA01FQT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O48500
Immunogen	Recombinant Escherichia phage T5 (Enterobacteria phage T5) DUT protein
Raised In	Rabbit
Species Reactivity	Escherichia phage T5 (Enterobacteria phage T5)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage T5 (Enterobacteria phage T5)
Lead Time	Made-to-order (14-16 weeks)
Target Names	DUT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.