



EEF1D Antibody

Product Code	CSB-PA007431XA01RB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53787
Immunogen	Recombinant Oryctolagus cuniculus (Rabbit) EEF1D protein
Raised In	Rabbit
Species Reactivity	Oryctolagus cuniculus (Rabbit)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Lead Time	Made-to-order (14-16 weeks)
Target Names	EEF1D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.